

Pocket Digital Alcohol Checker





PAL-Whiskey Cat.No. 4540

Reliable! Accuracy ±1.0%

PAL-Vodka

Cat.No. 4541

Fast! Only 3 seconds

PAL-Tequila

Cat.No. 4542

Small amount! Just 0.3ml

Lightl Only 100g




PAL- 焼酎 PAL-Gin Cat.No. 4543 Cat.No. 4544

PAL-Rum Cat.No. 4545



PAL-Brandy Cat.No. 4546



<Brewed beverage> Distill the beverage before measuring

Why can a refractometer measure the alcohol content of beverages?

A refractometer is a device that measures the total concentration of liquids using the refraction of light. As distilled beverages consist mainly of alcohol and water, the concentration of alcohol can be measured by a refractometer. Brewed beverages contain also sugar. Measuring them without distilling will give an inaccurate reading, as the value of sugar will be added to that of alcohol. Please distill brewed beverages before measuring.

Measurement Method





Apply 2-3 drops on the Press the START key.



Measurement value and prism temperature

Common Specifications

	Ethyl alcohol	
Measurement range	0.0 to 53.0%(ml/100ml)	
	5.0 to 30.0°C	
Resolution	0.5% / 0.1°C	
Accuracy	±1.0% / ±1°C	
Automatic temperature compensation range	5 to 30°C	

	Sample volume	At least 0.3ml
	Measurement time	Approx. 3 seconds
	Power supply	Two (2) AAA alkaline batteries
	Battery life	Approx. 11,000 (when using alkaline batteries)
	International Protection class	IP65
_	Dimensions and weight	55(W)×31(D)×109(H)mm, 100g (main unit only)

All ATAGO refractometers are designed and manufactured in Japan.

HACCP GMP GLP

ATAGO products comply with HACCF GMP, and GLP standards.



