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




ZEULAB
Making Testing Easier

Lab2Go

Set of tests for detection of
allergens in surfaces

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INTRODUCTION

LAB2GO is a set of rapid tests to detect traces of allergens on surfaces, presented in a briefcase that includes all the necessary components to perform the tests. This format facilitates the transport of the tests to the sampling site, either to perform the tests on site or to collect the sample and analyze it later.

They are immunochromatographic tests based on the detection of the antigen (allergen) that reacts in the strip with the specific antibody bound to red particles.

It offers personalized combinations between: gluten, egg, milk, soy, fish, crustacean, hazelnut and almond in a 10-test format per selected allergen.

COMPONENTS

STRIPS	10 / allergen
ANALYSIS BUFFER	10 / allergen
SWABS	10 / allergen
TEMPLATES (10x10 cm)	10 / allergen
PRODUCT CERTIFICATE	YES

After using all the components for each allergen, they can be separately purchased in order to complete the briefcase again. The references are shown in the table below:

Gluten (GL)	ZE/PR/GL10SW	Fish (P)	ZE/PR/P10SW
Egg (OV)	ZE/PR/OV10SW	Crustacean (CR)	ZE/PR/CR10SW
Milk (LBC)	ZE/PR/LBC10SW	Hazelnut (AV)	ZE/PR/AV10SW
Soy (S)	ZE/PR/S10SW	Almond (AL)	ZE/PR/AL10SW

STABILITY AND STORAGE

Store the test under dry conditions at room temperature (15-25°C). The kit has a shelf life of 12 months. Check the expiration date on the package.

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
NOTES

- The strips are very sensitive to moisture, keep them inside the supplied container and fully closed with the anti-humidity bag.
- Each allergen has its own buffer analysis associated with a specific colour. Please consult the legend provided in the last page.
- Use a swab, an analysis buffer tube and a strip for each sample.
- The test limit of detection is subjected to the use of the swabs provided.
- The milk test allows to detect the presence of milk with two proteins: β -lactoglobulin and casein. The double detection is particularly useful because only one test is required to identify the presence of the allergen, independently of its origin: milk, whey or caseinates.
- For manipulating the strip, hold it from the extreme opposite end of the arrows' area. It is recommended not to use gloves.
- The kit is designed for detecting allergen traces. The area to be analysed should pass a previous visual clean control.
- Select a representative area of at least 100 cm². The result obtained with the test indicates the presence or absence of the allergen on the analysed surface; it cannot be extrapolated into any value in ppm.
- Please, contact ZEULAB for further information.

PROCEDURE OF THE TEST


1. Open an ANALYSIS BUFFER tube, introduce a swab into it and soak the swab in the ANALYSIS BUFFER. The swab must be held and manipulated by its opposite end in order to avoid cross contamination.
2. Rub the swab through all possible directions on the analysed surface. Rotate from time to time the swab's point in contact with the surface using the provided template.
3. Place the swab in the initial tube with the ANALYSIS BUFFER and stir gently into the solution.
4. Discard the swab and introduce the strip test into the tube following the arrows' direction.
5. Wait for 10 min and read the results.

INTERPRETATION OF RESULTS

For a VALID TEST, the upper control line has to be visible after the incubation. When no control line is shown, the result must be considered INVALID. This may occur when the assay has not been performed properly or when the reagents are damaged. Repeat the assay with a new strip. Contact ZEULAB for further information.

NEGATIVE: Only the control line is observed in the white central zone of the strip test.

POSITIVE: In addition to control line, one or two test lines appear in the white central zone.

LAB2GO is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB does not assume any legal responsibility. To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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