

sIRoPAD




for fast analysis of plastic flakes, granules, pellets and crushed materials

With this fully automatic NIR-unit it is possible to identify up to 200 g non-carbonblack colored plastic parts directly as flakes, pellets, granulates or crushed material in an easy sampling way, being able to determine the quality, purity and the composition of the material very fast almost in ppm levels within ~15 min!



extended
InGaAs-



- * Plastics from household, packing and electric/electronics sections
- * Purity control of granulates, foils, pellets, flakes & grinded materials
- * Equivalent for PET-roasting test
- * Scanning time: full plate ~15min
- * Adjustable scanning area, plastic sizes smaller than 1 mm detectable
- * Scanning steps of 1-10 mm
- * Analysis and documentation
- * Pointing out the contamination

Technical Data:

- * Power Supply: 100-230 VAC, 50/60 Hz
- * 815 x 595 x 260 mm (height without NIR-camera)
- * 815 x 595 x 384 mm (height with NIR-camera)
- * Weight: 14 kg



ES France - Département Bio-tests & Industries
127 rue de Buzenval BP 26 - 92380 Garches



Tél. 01 47 95 99 90



e-mail : bio@es-france.com
Site Web : www.es-france.com


The sizes of the measuring area (templates) and the scanning parameters can be individually set depending on the material. Either as a single scan or as an automated batch with several scans with an individual area setting and parameters. Results can be saved manually or automatically.

siRoPAD scans the equally distributed samples on the ceramic reflection plate (max. area 50x33 cm) in an X/Y-pass way, measures and saves the analysis results with all information via USB-connection. The ICM-evaluation-software can be installed at any Windows-OS based computer. The software allows detailed spectra viewing, loading, saving and editing. This possibility helps to develop own analysis applications.



The analysis can be displayed graphically on the siRoPAD-touchscreen and external PC-display. Several measurement statistics can be done. With Zoom-functions, it is possible to do manual crosschecks in the map during and/or after the measurement at places under interest in detail.

All measurement files can be manually combined into analysis statistics which can be viewed in tabular or graphical form or exported as Excel-files. It enables quality controls over time periods, production processes or supplier qualities. Differences in the specific weight, thickness, etc. can also be corrected in the template settings.




Following polymer types are in the database:

PET * PE * PP * PVC * PS * ABS * PA * PC * PBT * SAN * PMMA * POM * PPO * PC+ABS * PVC+ABS * PC+PET * Cellulose * PE+PET * PE+PA * PP+PET * PLA * pet? * pvc?



Measurement results can be grouped, listed, arranged and printed out on the mini-plotter and/or at the connected PC-printer as an analysis report.



Individual After-scans can be carried out at any place. Questionable particles can be pointed out by a red-laser dot. This unique feature allows to review the detected contamination in detail.

