



G-DenPyc 3900

Gas Pycnometer Analyzer

Gold APP Instruments China

Lead You to Particle World Better



ES France - Département Bio-tests & Industries
127 rue de Buzenval BP 26 - 92380 Garches



Tél. 01 47 95 99 90
Fax. 01 47 01 16 22



e-mail : bio@es-france.com
Site Web : www.es-france.com

Gas Displacement Density Pycnometer

G-DenPyc 3900

Measurement for Foams, Powders, Solids, Films, Slurries, Coatings, Polymers, Calcination etc.

Working Principle

A gas pycnometer is a laboratory device used for measuring the density—or more accurately the volume—of solids, be they regularly shaped, porous or non-porous, monolithic, powdered, granular or in some way comminuted, employing some method of gas displacement and the volume. A gas pycnometer is also sometimes referred to as a helium pycnometer.

Gas expansion pycnometer is also known as constant volume gas pycnometer. The simplest type of gas pycnometer (due to its relative lack of moving parts) consists of two chambers, one (with a removable gas-tight lid) to hold the sample and a second chamber of fixed, known (via calibration) internal volume – referred to as the reference volume or added volume. The device additionally comprises a valve to admit a gas under pressure to one of the chambers, a pressure measuring device – usually a transducer – connected to the first chamber, a valved pathway connecting the two chambers, and a valved vent from the second of the chambers. In practice the sample may occupy either chamber that is gas pycnometers can be constructed such that the sample chamber is pressurized first, or such that it is the reference chamber that starts at the higher pressure.

Unmatchable Features


- ✚ Programmable logic controller (PLC) system obtains high integration and strong anti-interference, software operated fully automated analysis enables freely choose of multi experiment modes.
- ✚ Installed with a detachable filter in sample cell's bottom, can prevent samples be suctioned into manifolds; inputting gas from cavity bottom can avoid samples splash and software-controlled H-Sorb two-stage-stepping mode further ends splash happening.
- ✚ Provided standard 10ml sample cells, it is OK even have limited samples. Also can customize any volume sample cell from 5ml to 180ml. Exclusive G-DenPyc filling technology supports to choose different volumetric aluminum blocks to fill into cells instead changing other sizes cells according to samples' volume, this innovative method realizes minimal free space in cells and improve experiment accuracy.

- ✚ Patented V-Sorb monolithic manifolds system can improve sealing performance and reduce dead space largely, enhance system temperature's uniformity and anti-interference ability, all lead to a high accuracy and repeatability data.
- ✚ Programmable built-in Win10 system computer with a 13.3" display to real time monitor testing procedures, USB type external keyboard and mouse, RJ45 ethernet port is available for accessing to internet, also supports external display by RS232.
- ✚ The integrated design makes the whole machine compact and can be moved by one person. Meet any desktop environment.
- ✚ Hinged self-aligning lid, with easy and quick operation, can ensure high repeatable chamber volume in every testing.



A world of applications by G-DenPyc 3900

Gas pycnometer is used extensively for characterizing a wide variety of solids such as ceramic, catalysts, filter medium, nuclear fuel, oil & chemical industry, soil, fertilizer, carbon black, hard coke, fiber, minerals, pharmacy, cosmetics, cement, powdered foodstuff, desiccant (drying agent), powdered metal, ion exchange resin, silicon gel, alumina, titanium dioxide, solid foam etc.



- Analysis Method: gas displacement method (volume expansion method)
- Function: true density analysis
- Experimental Pressure: external vacuum pump is option, can adopt negative pressure (0-1Bar) or positive pressure (1-2Bar) two modes to analysis
- Accuracy: $\pm 0.02\%$ of reading, repeatability $\pm 0.01\%$, resolution can reach 0.0001g/ml
- Sample Ports: three samples analysis stations with provided 10ml standard sample tubes.
- Pressure Accuracy: imported high-precision pressure transducer, range from 0-3 Bar, accuracy can reach 0.04% F.S., stability 0.025% F.S.

- Adsorbate Gas: high purity He or N₂ (99.999%)
- Programmable embedded system with touchable display, equips USB external keyboard and mouse operation; the computer can be connected through RS232 communication mode; software can be run on an external computer for testing at the same time.

Net Dimensions

Length: 43cm
With: 35cm
Height: 37cm
Weight: 25kg



① 13.3" display with built-in PC

② 3 analysis ports



Software Advantages


G-DenPyc 3900 can be controlled by a built-in keypad or an external computer system. The whole analysis status can be continuously displayed in the 13.3 inch display.

Monoblock aluminum and thermostatic system can produce a minimal temperature floating ($\pm 0.1^\circ\text{C}$), help a lot for big mass samples whose density change dramatically with temperature.


G-DenPyc 3900 system not only support multiple analyses in a fixed pressure but also provide target pressure analysis within a select range.

Offered user selected reports:

- ✓ True density
- ✓ Summary report
- ✓ Sample log
- ✓ User defined tabular reports
- ✓ Equilibration report etc.




density report for ceramic sample



Data Reports

Equilibrium Pressure	Empty Cell Volume (ml)	Free Volume (ml)	True Density (g/ml)
51118.04	52.901100	37.744344	2.221099
47134.81	52.901100	37.800652	2.229382
45549.39	52.901100	37.802056	2.229589
46767.28	52.901100	37.791990	2.228103
50934.95	52.901100	37.791131	2.227977
47014.11	52.901100	37.798511	2.229065
51084.13	52.901100	37.786677	2.227320
47065.65	52.901100	37.793106	2.228268

data report tabular for battery material



true density distribution report



Headquarters

Gold APP Instruments Corp. Limited
Room 811, New Material Building, No. 7th,
Fenghui Mid Rd., Haidian Dist.,
Beijing 100094, P.R. China
Tel/WeChat/WhatsApp:
0086-18201085158
Email: sales@goldapp.com.cn
Web: www.goldapp.com.cn

Laboratory

Room 601, New Material Building, No. 7th, Fenghui Middle Rd.,
Haidian District, Beijing 100094,
P.R.China
Tel: +86-10-58711838
Fax: +86-10-58711838

Branch Offices

Gold APP Instruments (Nanjing) Corp. Limited
Room 512nd, No 4th Building,
Mingfa Commerce Square, No. 99th,
Yulan Rd., Yuhua District,
Nanjing 210012,
P.R.China
Tel: +86-25-58491095
Fax: +86-25-58491095

Gold iCON Instruments (Wuhan) Corp. Limited
Room 5068, No. 1st Building,
Huiyuan Block, No. 1st Rd.,
Wuhan University Science Park,
East Lake High-Tech Zone,
Wuhan 430223,
P.R.China
Tel: +86-27-59712850/1/2
Fax: +86-27-59712851 Ext.616

Our policy of continuous development may cause the information and specifications contained herein to change without notice or liability.

GOLD APP INSTRUMENTS®



ES France - Département Bio-tests & Industries
127 rue de Buzenval BP 26 - 92380 Garches



Tél. 01 47 95 99 90
Fax. 01 47 01 16 22



e-mail : bio@es-france.com
Site Web : www.es-france.com