

Quick A

HC-10
Thermal Conductivity Tester

EKO

www.eko-asia.com
www.eko-eu.com
www.eko-usa.com



Beyond Accuracy.

(ES) Equipements Scientifiques SA - Département Bio-tests & Industries - 127 rue de Buzenval BP 26 - 92380 Garches
Tél. 01 47 95 99 90 - Fax. 01 47 01 16 22 - e-mail: bio@es-france.com - Site Web: www.es-france.com

Quick Test & Easy Operation.

Measures within 60 seconds


Applications:

Quality assurance

Production

Material Qualification

Material research and development



Features


- Standalone portable Thermal Conductivity tester
- Suitable for homogeneous and VIP materials
- Simple and easy measurement by placing the sensor on the sample
- Very fast measurement within 60 seconds
- Capable of measuring a wide range of materials
- Software available enabling operation, calibration and analysis through a PC

Materials:

- VIP
- Rubber
- Plastics
- Ceramic
- Glass

Principle


In case of quick λ , the thermal conductivity is calculated by heat flow at the specimen surface layer for a short time non-static heating measurement. Heat loss is measured with the Sensor unit as Delta T which is inverse proportional to the material thermal conductivity.



Note: For measuring the thermal conductivity of a non-homogenous or layered sample a device which uses the the static method is recommended.

When the amount of heat loss is small, the sample surface temperature is high

When there is a large heat loss, the sample surface temperature is low




Measurement Procedure

Step1 :

Place the sensor unit on the cooling socket for 30 sec.


Step2 :

Place the sensor unit on the sample, then press the Start button to commence the measurment sequence. After 30 seconds the thermal conductivity will be displayed.



Evaluation of Measurement

For quick evaluation of a sample, the thermal conductivity threshold can be set. The measurement result will be evaluated as A,B or C and displayed on the LCD and LED indicators.



Specifications:

| | HC-10 | |
|---|---|---|
| Thermal Conductivity Range | VIP Sample: | 0.001 to 0.015 W/m•K |
| | Homogeneous sample: | 0.03 to 5.00 W/m•K |
| Method | Non-static transient Heat-up method | |
| Calibration | VIP Sample: | 3 or 4 samples of the same structure but different thermal conductivity |
| | Homogeneous sample: | 3 or 4 different standard samples |
| Repeatability | ± 5% | |
| Accuracy | Based on standard sample calibration | |
| Evaluation | Measurement of a sample can be classified (A,B or C) depending on the measurement result by setting a TC range (λ1& λ2) | |
| Display | Black and white LCD with back light | |
| Parameters | Thermal conductivity, Temperature, Sample class A,B,C | |
| Operating Temp. Range | +10 to +40 °C | |
| Data Storage | 99 measurement data and 20 calibration data | |
| Software | Data view and data management, Windows OS, English | |
| Interface | USB | |
| Standard sample (Attached to the Main Unit) | TEMPAX Glass, Acrylic and EPS | |
| Power, Power Consumption | AC adapter, AC100V to 240V 50/60Hz, DC24V Approx. 30W | |
| Size, Weight | W 250 x D 200 x H 85 mm, 4kg (Include Main Unit, Sensor Unit, and PowerSupply) | |

