

## **FLEXI ONE**

# **GPC Cleanup System**







#### Introduction .....

FLEXI ONE Gel Permeation Chromatography (GPC) is smart in size with modular design for easy-to-use concept. It could be applied to cleanup matrix such as food, grain, oil, herb, and other complicated samples from environment or biological tissues. Multiple standards like EPA (SW-846, 3640A), AOAC, S19 and Chinese national standards are applicable.

#### **Features**

- © Modular design with 8 inch touch screen with spectrogram displayed
- Controllable by laptop with 4 unites simultaneously
- Standard 250mL or 100mL receiving flask



### **Dedicated to Sample Preparation for 16 years**

#### Functions

- Embedded dual plunger pump with small pulse, low noise and broad velocity range
- Compatible with constant-wavelength UV detector or variable-wavelength UV detector
- Different column in size for various applications
- Quota-sampling by injection valve
- Automatic switch between sample components and effluent without watch-over
- ◆ Automatic caution for empty mobile phase solvent
- USB port for data export and printer connection
- Maximum 48 method store up
- Smart in size and favorable mobility



#### **Applications**



organo-chlorine pesticide



Phthalate



Pesticide residues & Veterinary drug residue



Dioxin

 $Statement: Specification\ may\ change\ without\ notice.\ Pree Kem\ reserves\ rights\ of\ final\ explanation.\ Version:\ 20160401$ 

Preekem Scientific instruments Co. Ltd.