

Sedimentation Device ZL 120

16

- The system is used for determining gluten quality and hazardous insects damage.
- LCD display with 6 lines.
- 6 timers
- 40 rpm induced motor enables system to perform the tests with standard 40 rpm rate, regardless of the load on the device.
- Every timer has a led being active when it is operated.
- The device enables operator to run only the timer without operating the motor.
- It stops automatically warning with a signal.
- It operates the test with standard angle.
- Possibility to measure both sedimentation and modified sedimentation value.



ERKAYA
INSTRUMENTS