## **Digital Viscometer RV Series**



## Feature

The viscometer is a new instrument designed from the circulating viscosity instrument, which has the advanced mechanical design and is controlled by microprocessor. It has the blue high light, blight LCD display can make you see clearly. The user can set the speed or spindle directly by press button and can get all the information about speed, spindle, temperature, viscosity from the disaply, including the percentage of the result in the whole range, aso the RV scan automatically select the spindle and speed by itself.



## Function

The viscometer has high sense, stable result, operating coveniently, beautiful outlook, which is used widely to measure the viscosity of Newton or Un-Newton liquid, such as adhesive solvent, lactop-rene, biolofical solvent, paint, oil, cosmetic, prinking ink and so on.

## Parameter

Test Range: RV-1S/1T (0.1-100,000 mpa.s) RV-2M/2T (0.1-2,000,000 mpa.s) Spindle: No. 0, 1, 2, 3, 4 Speed: RV-1S/1T (6, 12, 30, 60 rpm/min) RV-2M/2T (0.3, 0.6,1.5, 6, 12, 30, 60 rpm/min) RV-6M/6T (0.1/ 0.3 / 0.6 / 1.5 / 3 / 6 / 12 / 30 / 60 rpm/min) Accuracy: +/-1% Outline: 370\*325\*280 mm Net weight: 6.8kg

RV-6M/6T (0.1-6,000,000 mpa.s)

2



s 🛞

Tél. 01 47 95 99 90 Fax. 01 47 01 16 22