




SYM Noventa / 90°

for Microtome & Ultra Microtome

Easier and more free
section collection

Realize the work you want





Since the boat has a simple shape,
it is possible to use with your wafers and glass plates,
and set them up in various position.




Features of SYM Noventa

Point 1: Ingenuity in the boat shape




The boat section is the deepest compared to other Jumbo types. Then, the lower part of the blade is vertical and the corners are squared, so glass plates or wafers can be stably installed in the desired position.

Point 2: Design around the cutting edge



The area around the blade has a simple shape. This makes it possible to approach the wafers and glass plates close to the blade, making the collection easier.

Point 3: Ideal for SYM Cuesta and other manipulators



It is suitable for collection by hand, but it is more efficient when used in conjunction with any manipulators or SYM Cuesta.

The line-up of SYM Noventa

item No.	Width	Angle
NV2035	2.0mm	35°
NV2045	2.0mm	45°
NV3035	3.0mm	35°
NV3045	3.0mm	45°
NV4035	4.0mm	35°
NV4045	4.0mm	45°

More Freely and Easy to use !

