

Multi-Gen®

The better alternative to plastic disposable generator probes

Multi-Gen® generator probes are the better alternative to disposable plastic generator probes. This unique equipment is a reusable 7mm diameter generator probe manufactured from superior 316 stainless steel and PTFE, which allows for homogenizing samples without cross contamination and without the need to stop and clean a generator between samples. A Multi-Gen generator probe can quickly lock into the special Multi-Gen Adapter to provide hands-free ejection of the probe, minimizing the chance of accidental contamination. Multi-Gen generator probes are offered in multi-packs of 12 or 24 for convenience.

Advantages over plastic disposable generator probes:

First advantage, Multi-Gen generator probes can process more difficult samples, including tough tissues.

Second advantage, Multi-Gen generator probes can be sterilized by any method, including flaming and autoclaving.

Third advantage, Multi-Gen generator probes are more economical because of its ability to be reused an unlimited number of times not like plastic disposable generator probes

Multi-Gen 7 Features:

- Flat bottom design ideal for liquid and soft samples
- Shaft length best for microtubes to 5ml tubes

Multi-Gen 7XL Features:

- Open slotted tooth design ideal for the softest to the toughest of samples
- Longer shaft, approximately 115mm perfect for sample volumes from microtubes to 50ml tubes

Multi-Gen 7 and 7XL COMPONENTS AVAILABLE INDIVIDUALLY

ltem No.	Description
07-07200	Multi-Gen Adapter Required when using either Multi-Gen generator
02-070MG-12	Multi-Gen 7 Generators (12/package)
02-070MG-24	Multi-Gen 7 Generators (24/package)
02-070MGXL-12	Multi-Gen 7XL Generators (12/package)
80-00218	Multi-Gen Ejector Clamp Assembly For use with 80-00200 PR0200 Stand Assembly (see page 3) and Multi-Gen 7 generators only. Includes: Post clamp assembly, cross rod assembly, ejector collar
03-11703P	PTFE Bearings Replacement PTFE bearings for use in Multi- Gen 7 generators
Most recomme	nded motor unit for use with the Multi-Gen 7:
01-01200	PRO200 Homogenizer 120v homogenizer (see page 3)
01-02200	PRO200 Homogenizer 220v homogenizer (see page 3)
	PRO200 Stand Assembly 7 x 10 in. vinyl-padded hase with easy-life mechanism (see page 3)

NEXT-DAY DELIVERY ON THIS ITEM and most parts and accessories





