


# Roller Type Laboratory Mill RM 1300



- CE certified.
- Rolls are to grind the wheat, a specially designed sieve to separate the standard flour and the bran.
- 2 different collectors stocking flour and the bran separately.
- Performs the test with the flour efficiency of 55-70 % depending on the wheat variety and the moisture content.
- Possibility to grind dampened wheat up to 16-17 % moisture content.
- Possibility to see grinding by the glass front panel.
- Possibility to adjust the feeding rate according the moisture content or particle size of the wheat.
- 380 Volt 50 Hz. Electricity

